

# ABLIM1 Antibody



PACO07536

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ABLIM1

**Uniprot**

O14639

**Synonyms:**

actin binding LIM protein 1; ABLIM1; ABLIM; DKFZp781D0148; FLJ14564; KIAA0059; LIMAB1; LIMATIN; MGC1224 ;

**Immunogen:**

Human ABLIM1.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A