

ACAT1 Antibody



PACO07551

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:**Protein Background:****Gene ID:**

ACAT1

Uniprot

P24752

Synonyms:

acetyl-Coenzyme A acetyltransferase 1; ACAT1; ACAT; MAT; T2; THIL ;

Immunogen:

Human ACAT1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A