

ACE2 Antibody



PACO07556

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

ACE2

Uniprot

Q9BYF1

Synonyms:

angiotensin I converting enzyme (peptidyl-dipeptidase A) 2; ACE2; ACEH;
DKFZp434A014 ;

Immunogen:

Human ACE2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A