

ACLY Antibody



PACO07558

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

ACLY

Uniprot

P53396

Synonyms:

ATP citrate lyase; ACLY; ACL; ATPCL; CLATP ;

Immunogen:

Human ACLY.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A