

ACO2 Antibody



PACO07560

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ACO2

Uniprot

Q99798

Synonyms:

aconitase 2, mitochondrial; ACO2; ACONM; MGC20605; MGC33908 ;

Immunogen:

Human ACO2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A