

ACOT2 Antibody



PACO07562

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ACOT2

Uniprot

P49753

Synonyms:

acyl-CoA thioesterase 2; ACOT2; Mte1; PTE2; ZAP128 ;

Immunogen:

Human ACOT2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A