

ACOX2 Antibody



PACO07565

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:**Protein Background:****Gene ID:**

ACOX2

Uniprot

Q99424

Synonyms:

acyl-Coenzyme A oxidase 2, branched chain; ACOX2; BCOX; BRCACOX; BRCOX; THCCox ;

Immunogen:

Human ACOX2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A