

# ACP5 Antibody



PACO07567

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:****Protein Background:****Gene ID:**

ACP5

**Uniprot**

P13686

**Synonyms:**

acid, phosphatase 5, tartrate resistant; ACP5; MGC117378; TRAP ;

**Immunogen:**

Human ACP5.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A