

ACRV1 Antibody



PACO07570

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

ACRV1

Uniprot

P26436

Synonyms:

acrosomal vesicle protein 1; ACRV1; D11S4365; SP-10; SPACA2 ;

Immunogen:

Human ACRV1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A