

ACSL1 Antibody



PACO07571

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

ACSL1

Uniprot

P33121

Synonyms:

acyl-CoA synthetase long-chain family member 1; ACSL1; ACS1; FACL1; FACL2; LACS; LACS1; LACS2 ;

Immunogen:

Human ACSL1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A