

ACSS3 Antibody



PACO07576

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

ACSS3

Uniprot

Q9H6R3

Synonyms:

acyl-CoA synthetase short-chain family member 3; ACSS3; FLJ21963 ;

Immunogen:

Human ACSS3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A