

PACO07578

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat, Zebrafish

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

ACTB

Uniprot

P60709

Synonyms:

actin, beta; ACTB; PS1TP5BP1 ; beta-actin; -actin

Immunogen:

Human ACTB.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A