

ADCK5 Antibody



PACO07607

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ADCK5

Uniprot

Q3MIX3

Synonyms:

aarF domain containing kinase 5; ADCK5; FLJ35454; MGC126708 ;

Immunogen:

Human ADCK5.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A