

ADCY8 Antibody



PACO07611

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

ADCY8

Uniprot

P40145

Synonyms:

adenylate cyclase 8 (brain); ADCY8; AC8; ADCY3; HBAC1 ;

Immunogen:

Human ADCY8.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A