

ADRB2 Antibody



PACO07626

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

ADRB2

Uniprot

P07550

Synonyms:

adrenergic, beta-2-, receptor, surface; ADRB2; ADRB2R; ADRBR; B2AR; BAR; beta2AR ;

Immunogen:

Human ADRB2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A