

AGPAT9 Antibody



PACO07645

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

GPAT3

Uniprot

Q53EU6

Synonyms:

1-acylglycerol-3-phosphate O-acyltransferase 9; AGPAT9; AGPAT8; GPAT3; HMFN0839; LPAAT-theta; MAG1; MGC11324 ;

Immunogen:

Human AGPAT9.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A