

AGR3 Antibody



PACO07647

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

AGR3

Uniprot

Q8TD06

Synonyms:

anterior gradient homolog 3 (Xenopus laevis); AGR3; AG3; BCMP11; HAG3; PDIA18; hAG-3 ;

Immunogen:

Human AGR3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A