

# PHYKPL Antibody



PACO07651

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PHYKPL

**Uniprot**

Q8IUZ5

**Synonyms:**

alanine-glyoxylate aminotransferase 2-like 2; AGXT2L2; MGC117348; MGC15875; MGC45484 ;

**Immunogen:**

Human AGXT2L2.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A