

# ALDH2 Antibody



PACO07683

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:****Protein Background:****Gene ID:**

ALDH2

**Uniprot**

P05091

**Synonyms:**

aldehyde dehydrogenase 2 family (mitochondrial); ALDH2; ALDH-E2; ALDH1; ALDM; MGC1806 ;

**Immunogen:**

Human ALDH2.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A