

ALDOA Antibody



PACO07692

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

ALDOA

Uniprot

P04075

Synonyms:

aldolase A, fructose-bisphosphate; ALDOA; ALDA; MGC10942; MGC17716; MGC17767 ;

Immunogen:

Human ALDOA.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A