

ALKBH1 Antibody



PACO07699

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ALKBH1

Uniprot

Q13686

Synonyms:

alkB, alkylation repair homolog 1 (E. coli); ALKBH1; ABH; ABH1; ALKBH; alkB; hABH ;

Immunogen:

Human ALKBH1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A