

ALKBH5 Antibody



PACO07701

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ALKBH5

Uniprot

Q6P6C2

Synonyms:

alkB, alkylation repair homolog 5 (E. coli); ALKBH5; FLJ20308; OFOXD; OFOXD1 ;

Immunogen:

Human ALKBH5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A