

ALKBH7 Antibody



PACO07702

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ALKBH7

Uniprot

Q9BT30

Synonyms:

alkB, alkylation repair homolog 7 (E. coli); ALKBH7; MGC10974; SPATA11; UNQ6002 ;

Immunogen:

Human ALKBH7.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A