

AMOTL1 Antibody



PACO07719

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

AMOTL1

Uniprot

Q8IY63

Synonyms:

angiotensin like 1; AMOTL1; JEAP ;

Immunogen:

Human AMOTL1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C; C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A