

ANGPTL2 Antibody



PACO07732

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ANGPTL2

Uniprot

Q9UKU9

Synonyms:

angiopoietin-like 2; ANGPTL2; ARP2; HARP; MGC8889 ;

Immunogen:

Human ANGPTL2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A