

ANO8 Antibody



PACO07743

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

ANO8

Uniprot

Q9HCE9

Synonyms:

anoctamin 8; ANO8; KIAA1623; TMEM16H ;

Immunogen:

Human ANO8.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A