

ANP32A Antibody



PACO07744

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

ANP32A

Uniprot

P39687

Synonyms:

acid, c (leucine-rich) nuclear phosphoprotein 32 family, member A; ANP32A; C15orf1; I1PP2A; LANP; MAPM; MGC119787; MGC150373; PHAP1; PHAPI; PP32 ;

Immunogen:

Human ANP32A.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A