

ANXA5 Antibody



PACO07755

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ANXA5

Uniprot

P08758

Synonyms:

annexin A5; ANXA5; ANX5; ENX2; PP4 ;

Immunogen:

Human ANXA5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A