

AP3B1 Antibody



PACO07768

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

AP3B1

Uniprot

O00203

Synonyms:

adaptor-related protein complex 3, beta 1 subunit; AP3B1; ADTB3; ADTB3A; HPS; HPS2; PE ;

Immunogen:

Human AP3B1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A