

AP4M1 Antibody



PACO07773

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

AP4M1

Uniprot

O00189

Synonyms:

adaptor-related protein complex 4, mu 1 subunit; AP4M1; MU-4; MU-ARP2 ;

Immunogen:

Human AP4M1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A