

APBB2 Antibody



PACO07777

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

APBB2

Uniprot

Q92870

Synonyms:

amyloid beta (A4) precursor protein-binding, family B, member 2; APBB2; DKFZp434E033; FE65L; FE65L1; MGC35575 ;

Immunogen:

Human APBB2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A