

# APEH Antibody



PACO07779

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

APEH

**Uniprot**

P13798

**Synonyms:**

N-acylaminoacyl-peptide hydrolase; APEH; ACPH; APH; D3F15S2; D3S48E; DNF15S2; MGC2178; OPH ;

**Immunogen:**

Human APEH.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A