

APEX1 Antibody



PACO07780

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

APEX1

Uniprot

P27695

Synonyms:

APEX nuclease (multifunctional DNA repair enzyme) 1; APEX1; APE; APE1; APEN; APEX; APX; HAP1; REF1 ;

Immunogen:

Human APEX1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A