

APH1A Antibody



PACO07781

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

APH1A

Uniprot

Q96BI3

Synonyms:

anterior pharynx defective 1 homolog A (C. elegans); APH1A;6530402N02Rik; APH-1A; CGI-78 ;

Immunogen:

Human APH1A.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A