

NAA11 Antibody



PACO07804

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

NAA11

Uniprot

Q9BSU3

Synonyms:

ARD1 homolog B (S. cerevisiae); ARD1B; ARD2; MGC10646; hARD2 ;

Immunogen:

Human ARD1B.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A