

ARFGAP1 Antibody



PACO07810

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

ARFGAP1

Uniprot

Q8N6T3

Synonyms:

ADP-ribosylation factor GTPase activating protein 1; ARFGAP1; ARF1GAP; HRIHFB2281; MGC39924 ;

Immunogen:

Human ARFGAP1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A