

# ARHGAP10 Antibody



PACO07818

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ARHGAP10

**Uniprot**

A1A4S6

**Synonyms:**

Rho GTPase activating protein 10; ARHGAP10; FLJ20896; FLJ41791; GRAF2; PS-GAP; PSGAP ;

**Immunogen:**

Human ARHGAP10.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A