

ARHGAP9 Antibody



PACO07827

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ARHGAP9

Uniprot

Q9BRR9

Synonyms:

Rho GTPase activating protein 9; ARHGAP9;10C; FLJ16525; MGC1295; RGL1 ;

Immunogen:

Human ARHGAP9.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A