

ARHGEF10L Antibody



PACO07832

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ARHGEF10L

Uniprot

Q9HCE6

Synonyms:

Rho guanine nucleotide exchange factor (GEF) 10-like; ARHGEF10L; FLJ10521; GrinchGEF; KIAA1626; RP11-473A10.1 ;

Immunogen:

Human ARHGEF10L.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A