

ARL6IP5 Antibody



PACO07850

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ARL6IP5

Uniprot

O75915

Synonyms:

ADP-ribosylation-like factor 6 interacting protein 5; ARL6IP5; DERP11; GTRAP3-18; HSPC127; JWA; PRAF3; addicsin; hp22; jmx ;

Immunogen:

Human ARL6IP5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A