

# ARL8A Antibody



PACO07852

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

ARL8A

**Uniprot**

Q96BM9

**Synonyms:**

ADP-ribosylation factor-like 8A; ARL8A; ARL10B; FLJ45195; GIE2 ;

**Immunogen:**

Human ARL8A.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A