

ARMCX2 Antibody



PACO07857

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ARMCX2

Uniprot

Q7L311

Synonyms:

armadillo repeat containing, X-linked 2; ARMCX2; ALEX2; KIAA0512; MGC13343; MGC8742 ;

Immunogen:

Human ARMCX2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A