

# ARPC4 Antibody



PACO07862

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ARPC4

**Uniprot**

P59998

**Synonyms:**

actin related protein 2/3 complex, subunit 4, 20kDa; ARPC4; ARC20; MGC13544; p20-Arc ;

**Immunogen:**

Human ARPC4.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A