

ART3 Antibody



PACO07871

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ART3

Uniprot

Q13508

Synonyms:

ADP-ribosyltransferase 3; ART3; FLJ26404 ;

Immunogen:

Human ART3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A