

ART5 Antibody



PACO07872

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

ART5

Uniprot

Q96L15

Synonyms:

ADP-ribosyltransferase 5; ART5; MGC22848 ;

Immunogen:

Human ART5.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A