

ASB9 Antibody



PACO07876

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ASB9

Uniprot

Q96DX5

Synonyms:

ankyrin repeat and SOCS box-containing 9; ASB9; DKFZp564L0862; FLJ20636; MGC4954 ;

Immunogen:

Human ASB9.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A