

ASCC3 Antibody



PACO07879

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ASCC3

Uniprot

Q8N3C0

Synonyms:

activating signal cointegrator 1 complex subunit 3; ASCC3; ASC1p200; DJ467N11.1; HELIC1; MGC26074; RNAH; dJ121G13.4 ;

Immunogen:

Human ASCC3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A