

ASGR2 Antibody



PACO07883

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

ASGR2

Uniprot

P07307

Synonyms:

asialoglycoprotein receptor 2; ASGR2; ASGP-R; CLEC4H2; Hs.1259; L-H2 ;

Immunogen:

Human ASGR2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A