

PACO07889

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ASPH

**Uniprot**

Q12797

**Synonyms:**

aspartate beta-hydroxylase; ASPH; BAH; CASQ2BP1; HAAH; JCTN; junctin ;

**Immunogen:**

Human ASPH.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A