

ASZ1 Antibody



PACO07893

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ASZ1

Uniprot

Q8WWH4

Synonyms:

ankyrin repeat, SAM and basic leucine zipper domain containing 1; ASZ1; ALP1; ANKL1; C7orf7; GASZ; MGC26634; Orf3 ;

Immunogen:

Human ASZ1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A